

FIG. 1







FIG. 2





**FIG. 3A**

**FIG. 3B**






FIG. 3C






FIG. 3D





**FIG. 3E**



**FIG. 3F**






FIG. 3G




FIG. 3H






**Medtronic**

# MICROFRANCE<sup>®</sup>

## Josephson Suction Round Knife



The MicroFrance<sup>®</sup> Josephson Suction Round Knives were designed in conjunction with Gary D. Josephson, MD. By having suction on the knife plane, the line of sight stays clear of fluid without the surgeon having to hold a separate suction tube in the other hand. The surgeon can thus utilize more instruments simultaneously, allowing for quicker, more efficient otologic procedures.

For further information, please call Medtronic ENT at 800-874-5797 and/or consult Medtronic ENT's website at [www.MedtronicENT.com](http://www.MedtronicENT.com).

### ordering information

MCO5D	Josephson Suction Round Knife, 2mm	1 ea
MCO5E	Josephson Suction Round Knife, 3mm	1 ea

